



## CHANGELOG

# 3shape

VERSION	DATE OF CHANGE	DETAILS OF CHANGE
<b>v1.5</b>	14.12.2020	<p>Installation Guide separated and added to the zip-files of exo-cad, 3shape and Dental Wings, so they are better and easier reachable.</p> <p>Optimized all labanalog contours for 3D-printed jaw models. Now they all have a small Retentionstep, which de-forms when the labanalog is pushed into the printed orifice and holds the labanalog in its final position.</p> <p>Optimized the 3D-print contour for Safe-on-Four and Mini-Multi-Unit labanalog.</p> <p>Better interfaces for all bonding bases and Multi-Unit-Abutments for better visualization within the softwares.</p> <p>Assigned unique IDs for every part of the library so nothing will accidentally be overwritten by other libraries.</p> <p>Preface libraries with A-Line Milling Blanks for Arum milling machines have been postponed to the next update.</p>
<b>v1.4</b>	05.03.2020	<p>Added laboratory analog contours to the Multi-Unit-library</p>
<b>v1.3</b>	12.02.2020	<p>Added the laboratory analog contours to the Preface libraries</p>
<b>v1.2</b>	26.11.2019	<p>Added further compatibilities for the growing 3Shape Software Portfolio</p>
<b>v1.1</b>	19.08.2019	<p>Corrected scanbodies for SICvantage blue</p>
<b>v1.0</b>	04.03.2019	<p>Official Re-release</p>